

```

(***)Input identifying information(***)

projType = ToString[Evaluate[Input["Input the projection neuron type (i.e CC or CT)"]]];
mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]];
firstSection = ToString[
    Evaluate[Input["Input the number of the first histological section in the series"]]];
lastSection = ToString[
    Evaluate[Input["Input the number of the last histological section in the series"]]];
sectRange = Range[firstSection, lastSection];

(*****)

In[ ]:= laminarFluorPerSection = Table[Drop[Flatten[Import[
    StringJoin["S:/Imaging/Garrett/BCMM_EpiFluorConfocScope/Other/V1_CCandCTcells/",
    projType, "/", mouse, "/LaminarFluorData/", mouse, "_Section",
    ToString[sectRange[[n]]], ".xlsx"], 1], 2], {n, 1, Length[sectRange]}];

In[ ]:= resol = (Last[laminarFluorPerSection[[1]][[All, 1]]] -
    First[laminarFluorPerSection[[1]][[All, 1]]]) /
    Length[laminarFluorPerSection[[1]][[All, 1]]];

In[ ]:= (***)Median-smooth at 5 um(***)

In[ ]:= window = Round[5/resol];

In[ ]:= medianSmoothedPerSection =
    Table[Partition[Riffle[(laminarFluorPerSection[[n]])[[All, 1]],
        MovingMap[Quantile[#, 1/2] &, (laminarFluorPerSection[[n]])[[All, 2]],
        window, "Reflected"]], 2], {n, 1, Length[sectRange]}];

(***)Check that all sections go up to at least 150 um(***)

In[ ]:= Manipulate[ListLinePlot[medianSmoothedPerSection[[n]], PlotRange -> All],
    {n, 1, Length[sectRange], 1}]

In[ ]:= Table[Last[medianSmoothedPerSection[[n]][[All, 1]]], {n, 1, Length[sectRange]}]

In[ ]:= (***)Truncate all lists at 150 um(***)

In[ ]:= truncPoint = Round[150/resol];

In[ ]:= medianSmoothedPerSectionTrunc =
    Table[Take[medianSmoothedPerSection[[n]], truncPoint], {n, 1, Length[sectRange]}];

In[ ]:= (***)Normalize all lists to max fluorescence(***)

In[ ]:= maxFluor = Max[Flatten[
    Table[medianSmoothedPerSectionTrunc[[n]][[All, 2]], {n, 1, Length[sectRange]}]]];

In[ ]:= medianSmoothedPerSectionTruncNorm =
    Table[Partition[Riffle[medianSmoothedPerSectionTrunc[[n]][[All, 1]],
        medianSmoothedPerSectionTrunc[[n]][[All, 2]]/maxFluor],
    2], {n, 1, Length[sectRange]}];

```

```

In[ ]:= (**Calculate the mean laminar fluorescence profile for this animal and export data**)

In[ ]:= meanFluorProfile =
    Mean[Table[medianSmoothedPerSectionTruncNorm[[n]], {n, 1, Length[sectRange]}]];

    (**Export layer 1 profile of the mean
    normalized dendritic fluorescence from this animal**)

In[ ]:= Export[StringJoin["S:/Imaging/Garrett/BCMM_EpiFluorConfocScope/Other/V1_CCandCTcells/",
    projType, "/", mouse, "/", mouse, "_meanLaminarFluor.txt"], meanFluorProfile];

```