```
In[@]:= (***Input identifying information***)
In[=]:= date = ToString[Evaluate[Input["Input the date of the experiment"]]]
In[*]:= mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]]
Info |= sessionNumBef = Evaluate[Input["Input the session number before manipulation"]]
ln[\cdot]:= sessionNumAft = Evaluate[Input["Input the session number after manipulation"]]
In[*]:= discROIsQ = ToString[Evaluate[Input["Are there any non-useable ROIs?"]]]
In[@]:= numROIs =
       Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
           mouse, "/Session", ToString[sessionNumBef], "/dFOverF0TimeSeries/"]]]];
In[*]:= If[discROIsQ == "Yes",
     nonUseROIs = ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208 2PRig/",
           date, "/", mouse, "/Session", ToString[sessionNumAft], "/", date, " ", mouse,
           "_Session", ToString[sessionNumAft], "_nonUseableROIs.txt"]]];, nonUseROIs = {};]
Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
         mouse, "/Session", ToString[sessionNumBef], "/LocomotionData/", date, "_
         mouse, "_Session", ToString[sessionNumBef], "_locModROIs.txt"], "List"];
ln[e]:= sigROIsAft = Import[StringJoin["S:/Imaging/Garrett/FMB208 2PRig/", date, "/",
         mouse, "/Session", ToString[sessionNumAft], "/LocomotionData/", date, "_"
         mouse, "_Session", ToString[sessionNumAft], "_locModROIs.txt"], "List"];
Info |:= sigRespROIs = DeleteCases [Table [
         If[MemberQ[sigROIsBef, n] && MemberQ[sigROIsAft, n], n, Null], {n, 1, numROIs}], Null];
      (***ROIs are only acceptable if they were acceptable both before and after***)
Info |:= nonSigRespROIs = Complement[Range[numROIs], sigRespROIs]
In[e]= noGoodROIs = DeleteDuplicates[Join[nonSigRespROIs, nonUseROIs]]
ln[*]:= usefulROIs = Complement[Range[numROIs], noGoodROIs]
<code>ln[*]:= (***For each ROI that can be used for paired analysis,</code>
    upload the locomotion modulation index for the 2 sessions,
    and pair the indices for the sessions**)
m[\cdot]:= Table[Evaluate@ToExpression[StringJoin["locModBef", ToString[n]]] = ToExpression/@
          Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse, "/Session",
            ToString[sessionNumBef], "/LocomotionData/", date, "_", mouse, "_Session",
            ToString[sessionNumBef], "_", "SummaryLocModIndex_Baseline15sAway_ROI",
            ToString[n], ".txt"], "List"];, {n, usefulROIs}];
ln[*]:= Table[Evaluate@ToExpression[StringJoin["locModAft", ToString[n]]] = ToExpression/@
          Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse, "/Session",
            ToString[sessionNumAft], "/LocomotionData/", date, "_", mouse, "_Session",
            ToString[sessionNumAft], "_", "SummaryLocModIndex_Baseline15sAway_ROI",
            ToString[n], ".txt"], "List"];, {n, usefulROIs}];
```

```
Im[@]= Table[Evaluate@ToExpression[StringJoin["locModPaired", ToString[n]]] =
               {n, Flatten@{ToExpression[StringJoin["locModBef", ToString[n]]],
                    ToExpression[StringJoin["locModAft", ToString[n]]]}};, {n, usefulROIs}];
<code>ln[*]:= (***For each ROI that can be used for paired analysis,</code>
       upload the peri-loc-onset Z-scored dF/F for the 2 sessions**)
ln[*]:= Table[Evaluate@ToExpression[StringJoin["periOnTraceBef", ToString[n]]] =
              ToExpression /@ Import [StringJoin ["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
                    mouse, "/Session", ToString[sessionNumBef], "/LocomotionData/", date, "_", mouse,
                    "_Session", ToString[sessionNumBef], "_", "PeriOnsetZDFF_PreAndPostBaseline_ROI",
                    ToString[n], ".txt"], "List"];, {n, usefulROIs}];
log_{log} = Table[Evaluate@ToExpression[StringJoin["periOnTraceAft", ToString[n]]] = log_{log} = log
              ToExpression /@ Import [StringJoin ["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
                    mouse, "/Session", ToString[sessionNumAft], "/LocomotionData/", date, " ", mouse,
                    "_Session", ToString[sessionNumAft], "_", "PeriOnsetZDFF_PreAndPostBaseline_ROI",
                    ToString[n], ".txt"], "List"];, {n, usefulROIs}];
ln[∘]:= (***For each ROI that can be used for paired analysis,
       upload the peri-loc-offset Z-scored dF/F for the 2 sessions**)
In[*]:= Table[Evaluate@ToExpression[StringJoin["periOffTraceBef", ToString[n]]] =
              ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
                    mouse, "/Session", ToString[sessionNumBef], "/LocomotionData/", date, "_", mouse,
                    "_Session", ToString[sessionNumBef], "_", "PeriOffsetZDFF_PreAndPostBaseline_ROI",
                    ToString[n], ".txt"], "List"];, {n, usefulROIs}];
ln[@]:= Table[Evaluate@ToExpression[StringJoin["periOffTraceAft", ToString[n]]] =
              ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
                    mouse, "/Session", ToString[sessionNumAft], "/LocomotionData/", date, "_", mouse,
                    "_Session", ToString[sessionNumAft], "_", "PeriOffsetZDFF_PreAndPostBaseline_ROI",
                    ToString[n], ".txt"], "List"];, {n, usefulROIs}];
ln[*]: Manipulate[{ListLinePlot[{ToExpression[StringJoin["periOnTraceBef", ToString[n]]],
              ToExpression[StringJoin["periOnTraceAft", ToString[n]]]},
             PlotStyle → {Black, Green}, PlotRange → All],
           ListLinePlot[{ToExpression[StringJoin["periOffTraceBef", ToString[n]]],
               ToExpression[StringJoin["periOffTraceAft", ToString[n]]]},
             PlotStyle → {Black, Green}, PlotRange → All],
           ToExpression[StringJoin["locModPaired", ToString[n]]]}, {n, usefulROIs}]
ln[*]: Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
             "/PairedAnalysis/", date, "_", mouse, "_pairedROIsLoc.txt"], usefulROIs];
"/PairedAnalysis/", date, "_", mouse, "_locModPaired_ROI", ToString[n], ".txt"],
             ToExpression[StringJoin["locModPaired", ToString[n]]]], {n, usefulROIs}];
```