

```

In[ ]:= (**Input identifying information**)

In[ ]:= date = ToString[Evaluate[Input["Input the date of the experiment"]]]

In[ ]:= mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]]

In[ ]:= sessionNumBef = Evaluate[Input["Input the session number before manipulation"]]

In[ ]:= sessionNumAft = Evaluate[Input["Input the session number after manipulation"]]

In[ ]:= discROIsQ = ToString[Evaluate[Input["Are there any non-useable ROIs?"]]]

In[ ]:= numROIs =
    Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
        mouse, "/Session", ToString[sessionNumBef], "/dFOverF0TimeSeries/"]]]];

In[ ]:= If[discROIsQ == "Yes",
    nonUserROIs = ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/",
        date, "/", mouse, "/Session", ToString[sessionNumAft], "/", date, "_", mouse,
        "_Session", ToString[sessionNumAft], "_nonUseableROIs.txt"]]]]; nonUserROIs = {};

In[ ]:= sigROIsBef = Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
    "/Session", ToString[sessionNumBef], "/VisStimResults/", date, "_", mouse,
    "_Session", ToString[sessionNumBef], "_sigResponsiveROIs.txt"], "List"];

In[ ]:= sigROIsAft = Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
    "/Session", ToString[sessionNumAft], "/VisStimResults/", date, "_", mouse,
    "_Session", ToString[sessionNumAft], "_sigResponsiveROIs.txt"], "List"];

In[ ]:= sigRespROIs = DeleteCases[Table[
    If[MemberQ[sigROIsBef, n] || MemberQ[sigROIsAft, n], n, Null], {n, 1, numROIs}], Null];
(**An ROI is significantly responsive it was significantly responsive
during the before session OR the after session**)

In[ ]:= nonSigRespROIs = Complement[Range[numROIs], sigRespROIs]

In[ ]:= noGoodROIs = DeleteDuplicates[Join[nonSigRespROIs, nonUserROIs]]

In[ ]:= usefulROIs = Complement[Range[numROIs], noGoodROIs]

In[ ]:= (*****For each ROI picked for the session,
upload the average evoked dF/F values*****)

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFAvgBef", ToString[n]]] =
    ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date,
        "/", mouse, "/Session", ToString[sessionNumBef], "/VisStimResults/",
        date, "_", mouse, "_Session", ToString[sessionNumBef], "_",
        "overallVisDFFZScored_ROI", ToString[n], ".txt"], "List"]]; {n, 1, numROIs}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFAvgAft", ToString[n]]] =
    ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date,
        "/", mouse, "/Session", ToString[sessionNumAft], "/VisStimResults/",
        date, "_", mouse, "_Session", ToString[sessionNumAft], "_",
        "overallVisDFFZScored_ROI", ToString[n], ".txt"], "List"]]; {n, 1, numROIs}];

In[ ]:= (**Create a paired before-after data list for the evoked dF/F**)

```

```

In[ ]:= meanDFFpairs = Table[{(ToExpression[StringJoin["dFFAvgBef", ToString[n]])[[1]],
  (ToExpression[StringJoin["dFFAvgAft", ToString[n]])[[1]]}, {n, usefulROIs}];

```

```

In[ ]:= (*****)

```

```

In[ ]:= Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
  "/PairedAnalysis/", date, "_", mouse, "_visRespPaired.txt"], meanDFFpairs];

```