

```

(***)Input identifying information(***)

projType = ToString[
  Evaluate[Input["Input the projection type (i.e V1toPM, LMtoPM, or LPtoPM)"]]];

mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]];

firstSection = ToString[
  Evaluate[Input["Input the number of the first histological section in the series"]]];

lastSection = ToString[
  Evaluate[Input["Input the number of the last histological section in the series"]]];

sectRange = Range[firstSection, lastSection];

(*****)

In[ ]:= laminarFluorPerSection = Table[Join[Drop[Flatten[
  Import[StringJoin["S:/Imaging/Garrett/BCMM_EpifluorConfocScope/Other/
    LaminarAxonFluorescenceHistology/",
    projType, "/", mouse, "/LaminarFluorData/", mouse, "_Section",
    ToString[sectRange[[n]]], "_Upper.xlsx"], 1], 2],
  Partition[Riffle[(Drop[Flatten[Import[StringJoin[
    "S:/Imaging/Garrett/BCMM_EpifluorConfocScope/Other/
    LaminarAxonFluorescenceHistology/",
    projType, "/", mouse, "/LaminarFluorData/", mouse,
    "_Section", ToString[sectRange[[n]]], "_Lower.xlsx"], 1],
  2][[All, 1]] + Last[Drop[Flatten[Import[StringJoin[
    "S:/Imaging/Garrett/BCMM_EpifluorConfocScope/Other/
    LaminarAxonFluorescenceHistology/",
    projType, "/", mouse, "/LaminarFluorData/", mouse, "_Section",
    ToString[sectRange[[n]]], "_Upper.xlsx"], 1],
  2][[All, 1]]], (Drop[Flatten[Import[StringJoin[
    "S:/Imaging/Garrett/BCMM_EpifluorConfocScope/Other/
    LaminarAxonFluorescenceHistology/",
    projType, "/", mouse, "/LaminarFluorData/", mouse,
    "_Section", ToString[sectRange[[n]]], "_Lower.xlsx"], 1],
  2][[All, 2]]], 2]], {n, 1, Length[sectRange]}]];

In[ ]:= resol = (Last[laminarFluorPerSection[[1]][[All, 1]]] -
  First[laminarFluorPerSection[[1]][[All, 1]]]) /
  Length[laminarFluorPerSection[[1]][[All, 1]]];

In[ ]:= (***)Median-smooth at 10 um(***)

In[ ]:= window = Round[10/resol];

In[ ]:= medianSmoothedPerSection =
  Table[Partition[Riffle[(laminarFluorPerSection[[n]])[[All, 1]],
    MovingMap[Quantile[#, 1/2] &, (laminarFluorPerSection[[n]])[[All, 2]],
    window, "Reflected"]], 2], {n, 1, Length[sectRange]}];

(***)Check that all sections go up to at least 800 um(***)

```

```

In[ ]:= Manipulate[ListLinePlot[medianSmoothedPerSection[[n]], PlotRange → All],
  {n, 1, Length[sectRange], 1}]

In[ ]:= Table[Last[medianSmoothedPerSection[[n]][[All, 1]]], {n, 1, Length[sectRange]}]

In[ ]:= (***Truncate all lists at 800 um***)

In[ ]:= truncPoint = Round[800/resol];

In[ ]:= medianSmoothedPerSectionTrunc =
  Table[Take[medianSmoothedPerSection[[n]], truncPoint], {n, 1, Length[sectRange]}];

In[ ]:= (***Normalize all lists to max fluorescence***)

In[ ]:= maxFluor = Max[Flatten[
  Table[medianSmoothedPerSectionTrunc[[n]][[All, 2]], {n, 1, Length[sectRange]}]]];

In[ ]:= medianSmoothedPerSectionTruncNorm =
  Table[Partition[Riffle[medianSmoothedPerSectionTrunc[[n]][[All, 1]],
    medianSmoothedPerSectionTrunc[[n]][[All, 2]]/maxFluor],
    2], {n, 1, Length[sectRange]}];

In[ ]:= (***Calculate the mean laminar fluorescence profile for this animal and export data***)

In[ ]:= meanFluorProfile =
  Mean[Table[medianSmoothedPerSectionTruncNorm[[n]], {n, 1, Length[sectRange]}]];

  (***)Export laminar profile of the mean
  normalized axonal fluorescence from this animal***)

In[ ]:= Export[StringJoin[
  "S:/Imaging/Garrett/BCMM_EpifluorConfocScope/Other/LaminarAxonFluorescenceHistology/"
  , projType, "/", mouse, "/", mouse, "_meanLaminarFluor.txt"],
  meanFluorProfile];

```