

```

    (**Input identifying information**)

In[ ]:= date = ToString[Evaluate[Input["Input the date of the experiment"]]]
Out[ ]:= 071223

In[ ]:= mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]]
Out[ ]:= Mouse23666

In[ ]:= sessionNum = Evaluate[Input["Input the session number"]]
Out[ ]:= 2

In[ ]:= (**Import pupil diameter time series**)

In[ ]:= pupilTimeSeriesBef =
    ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
        mouse, "/Session", ToString[sessionNum], "/Pupil/", date, "_", mouse, "_",
        "Session", ToString[sessionNum], "_PupilRadiusTimeSeries.txt"], "List"]];

In[ ]:= timeVals = Part[#, 1] & /@ pupilTimeSeriesBef;

In[ ]:= (**Create pupil diameter time series,
    excluding NaN elements from edge detection abnormalities**)

In[ ]:= pupilTimeSeries =
    DeleteCases[Table[If[NumberQ@pupilTimeSeriesBef[[n, 2]], {pupilTimeSeriesBef[[n, 1]],
        pupilTimeSeriesBef[[n, 2]]}, Null], {n, 1, Length[pupilTimeSeriesBef]}], Null];

In[ ]:= (**Interpolate the pupil diameter time series and resample at 10 Hz**)

In[ ]:= pupilInterp = Interpolation[pupilTimeSeries];

In[ ]:= pupil10Hz = Table[pupilInterp[t],
    {t, First[pupilTimeSeries][[1]], Last[pupilTimeSeries][[1]], 0.1}];

In[ ]:= pupil10HzTimeSeries = Partition[Riffle[Table[n,
    {n, First[pupilTimeSeries][[1]], Last[pupilTimeSeries][[1]], 0.1}], pupil10Hz], 2];

In[ ]:= (**Further,
    only keep the dilations and constrictions that occur during non-locomotion,
    non-whisking periods**)

In[ ]:= quiescentPeriods =
    ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
        mouse, "/Session", ToString[sessionNum], "/LocomotionData/", date, "_", mouse,
        "_", "Session", ToString[sessionNum], "_QuiescentBouts.txt"], "List"]];

In[ ]:= (**Import calcium fluorescence traces**)

In[ ]:= numAdditionalMovies =
    Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/",
        date, "/", mouse, "/Session", ToString[sessionNum], "/Ftraces/"]]]] - 2;
    (*Subtract 2 because the first movie was already imported and the
    other file in the directory is an ROI list*)

```

```

In[ ]:= tpFrameTimes =
  Drop[Drop[(Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
    mouse, "/Session", ToString[sessionNum], "/", date, "_", mouse, "_",
    "Session", ToString[sessionNum], "_2PFrameTimes.txt"], "List"]], 16], -1];

In[ ]:= numFramesPerAcq = Round[Length[tpFrameTimes] / (numAdditionalMovies + 1)];

In[ ]:= tpFrameRate = Round[Mean[numFramesPerAcq / Flatten[Differences /@
  ({First[#], Last[#]} & /@ (Partition[tpFrameTimes, numFramesPerAcq]))]]];

In[ ]:= (**For each ROI picked for the session, upload the extracted dF/F0 time series**)

In[ ]:= numROIs =
  Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
    mouse, "/Session", ToString[sessionNum], "/dFoverF0TimeSeries/"]]]];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFtimeseries", ToString[n]]] =
  ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date,
    "/", mouse, "/Session", ToString[sessionNum], "/dFoverF0TimeSeries/",
    date, "_", mouse, "_Session", ToString[sessionNum], "_",
    "dFoverF0ts_ROI", ToString[n], ".txt"], "List"]; {n, 1, numROIs}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFts", ToString[n]]] =
  Partition[Riffle[(ToExpression[StringJoin["dFFtimeseries", ToString[n]]])][[All, 1]],
    Standardize[(ToExpression[StringJoin["dFFtimeseries", ToString[n]]])][[All, 2]]],
    2]; {n, 1, numROIs}];

  (*****
  *****)
  (**In this part, calculate cross-correlograms between dF/F0 and pupil diameter**)
  (*****
  *****)

In[ ]:= (**Interpolate the 10 Hz sampled pupil time series**)

In[ ]:= pupil10HzInt = Interpolation[pupil10HzTimeSeries];

In[ ]:= (**Lowpass filter the dF/F0 traces at 10 Hz**)

In[ ]:= Table[
  Evaluate@ToExpression[StringJoin["dFFlpfTS", ToString[n]]] = LowpassFilter[TimeSeries@
    ToExpression[StringJoin["dFFts", ToString[n]]], 10 Hz]; {n, 1, numROIs}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFlpfInt", ToString[n]]] =
  Interpolation[ToExpression[StringJoin["dFFlpfTS", ToString[n]]]; {n, 1, numROIs}];

In[ ]:= (**Only keep quiescent periods at least 9 s long**)

In[ ]:= quiescentPeriodsLongBef =
  DeleteCases[Table[If[quiescentPeriods[[n, 2]] - quiescentPeriods[[n, 1]] ≥ 9,
    quiescentPeriods[[n]], Null], {n, 1, Length[quiescentPeriods]}], Null];

In[ ]:= (**For each acquisition bout, determine the quiescent periods
  contained within it and then make sure they're at least 9 s long**)

```

```

In[ ]:= acqBoutIntervals =
  Interval /@ ({First[#] + 8, Last[#] - 8} & /@ (Partition[tpFrameTimes, numFramesPerAcq]));

In[ ]:= testedIntervals = Table[MinMax /@ (IntervalIntersection[acqBoutIntervals[[n]], #] & /@
  (Interval /@ quiescentPeriodsLongBef)), {n, 1, Length[acqBoutIntervals]}];

In[ ]:= acceptedQuiescentPeriods = Flatten[
  Table[DeleteCases[Table[If[NumberQ@ (Differences[testedIntervals[[m, n]]][[1]]),
    testedIntervals[[m, n]], Null], {n, 1, Length[testedIntervals[[m]]}],
    Null], {m, 1, Length[testedIntervals]}], 1];

In[ ]:= quiescentPeriodsLong = DeleteCases[
  Table[If[acceptedQuiescentPeriods[[n, 2]] - acceptedQuiescentPeriods[[n, 1]] ≥ 9,
    acceptedQuiescentPeriods[[n]], Null],
    {n, 1, Length[acceptedQuiescentPeriods]}], Null];

In[ ]:= (**Compute cross-correlations between pupil diameter and dF/F0 during quiescent
  periods. Compute in a time window of 8 seconds with 100 ms resolution**)

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFpupilCrossCorr", ToString[roi]]] =
  Mean[DeleteCases[Table[pupil = Table[pupil10HzInt[t],
    {t, quiescentPeriodsLong[[i, 1]], quiescentPeriodsLong[[i, 2]], 0.1}];
    neuron = Table[(ToExpression[StringJoin["dFFlpfInt", ToString[roi]]])[t],
    {t, quiescentPeriodsLong[[i, 1]], quiescentPeriodsLong[[i, 2]], 0.1}];
    If[quiescentPeriodsLong[[i, 2]] - quiescentPeriodsLong[[i, 1]] > 8,
      Join[Reverse[Table[Correlation[pupil, RotateRight[neuron, -n]], {n, 0, 80, 1}]],
        Drop[Table[Correlation[pupil, RotateRight[neuron, n]], {n, 0, 80, 1}], 1], Null],
    {i, 1, Length[quiescentPeriodsLong]}], Null];, {roi, 1, numROIs}];

In[ ]:= crossCorrTimeVals = Table[n, {n, -8, 8, 0.1}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFpupilCrossCorrTS", ToString[n]]] =
  Partition[Riffle[crossCorrTimeVals, (ToExpression[
    StringJoin["dFFpupilCrossCorr", ToString[n]])], 2];, {n, 1, numROIs}];

(**Visualize pupil-dFF cross correlograms**)

In[ ]:= Manipulate[ListLinePlot[ToExpression[StringJoin["dFFpupilCrossCorrTS", ToString[roi]]],
  {roi, 1, numROIs, 1}]

In[ ]:= (**Export dFF-pupil cross-correlation data**)

In[ ]:= Table[Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
  "/Session", ToString[sessionNum], "/", "Pupil/", date, "_", mouse, "_",
  "Session", ToString[sessionNum], "_dFFpupilCrossCorr_ROI", ToString[n], ".txt"],
  ToExpression[StringJoin["dFFpupilCrossCorrTS", ToString[n]]], {n, 1, numROIs}];

In[ ]:= (*****)

pupilModROIs = Range[numROIs];

In[ ]:= Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
  "/Session", ToString[sessionNum], "/", "Pupil/", date, "_", mouse, "_",
  "Session", ToString[sessionNum], "_pupilModROIs", ".txt"], pupilModROIs];

```