```
In[@]:= (***Input identifying information***)
In[=]:= date = ToString[Evaluate[Input["Input the date of the experiment"]]]
In[*]:= mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]]
log_{ij} = sessionNumBef = Evaluate[Input["Input the session number before manipulation"]]
ln[\cdot]:= sessionNumAft = Evaluate[Input["Input the session number after manipulation"]]
In[*]:= discROIsQ = ToString[Evaluate[Input["Are there any non-useable ROIs?"]]]
In[*]:= numROIs =
       Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
           mouse, "/Session", ToString[sessionNumBef], "/dFOverF0TimeSeries/"]]]];
In[*]:= If[discROIsQ == "Yes",
     nonUseROIs = ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/",
           date, "/", mouse, "/Session", ToString[sessionNumAft], "/", date, " ", mouse,
           "_Session", ToString[sessionNumAft], "_nonUseableROIs.txt"]]];, nonUseROIs = {};]
ln[*]: sigROIsBef = Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date,
         "/", mouse, "/Session", ToString[sessionNumBef], "/Pupil/", date, "_",
         mouse, "_Session", ToString[sessionNumBef], "_pupilModROIs.txt"], "List"];
ln[*]: sigROIsAft = Import[StringJoin["S:/Imaging/Garrett/FMB208 2PRig/", date,
         "/", mouse, "/Session", ToString[sessionNumAft], "/Pupil/", date, "_",
         mouse, "_Session", ToString[sessionNumAft], "_pupilModROIs.txt"], "List"];
In[@]:= sigRespROIs = DeleteCases[Table[
         If[MemberQ[sigROIsBef, n] && MemberQ[sigROIsAft, n], n, Null], {n, 1, numROIs}], Null];
      (***ROIs are only acceptable if they were acceptable both before and after***)
In[*]:= nonSigRespROIs = Complement[Range[numROIs], sigRespROIs]
In[*]: noGoodROIs = DeleteDuplicates[Join[nonSigRespROIs, nonUseROIs]]
ln[*]:= usefulROIs = Complement[Range[numROIs], noGoodROIs]
In[*]:= (***For each ROI that can be used for paired analysis,
    upload the locomotion modulation index for the 2 sessions,
    and pair the indices for the sessions**)
ln[*]:= Table[Evaluate@ToExpression[StringJoin["dffPupCCBef", ToString[n]]] =
         ToExpression /@ Import [StringJoin ["S:/Imaging/Garrett/FMB208_2PRig/",
            date, "/", mouse, "/Session", ToString[sessionNumBef], "/Pupil/",
            date, "_", mouse, "_Session", ToString[sessionNumBef], "_",
             "dFFpupilCrossCorr_ROI", ToString[n], ".txt"], "List"];, {n, usefulROIs}];
ln[@]:= Table[Evaluate@ToExpression[StringJoin["dffPupCCAft", ToString[n]]] =
         ToExpression /@ Import [StringJoin ["S:/Imaging/Garrett/FMB208_2PRig/",
            date, "/", mouse, "/Session", ToString[sessionNumAft], "/Pupil/",
            date, "_", mouse, "_Session", ToString[sessionNumAft], "_",
            "dFFpupilCrossCorr_ROI", ToString[n], ".txt"], "List"];, {n, usefulROIs}];
Info ]:= catenatedPupilCrossCorrsBef =
       Table[ToExpression[StringJoin["dffPupCCBef", ToString[n]]], {n, usefulROIs}];
```

```
Im[=::= peakCrossCorrsBef = Table[Max[(catenatedPupilCrossCorrsBef[[n]])[[All, 2]]],
       {n, 1, Length[catenatedPupilCrossCorrsBef]}];
Infol:= catenatedPupilCrossCorrsAft =
      Table[ToExpression[StringJoin["dffPupCCAft", ToString[n]]], {n, usefulROIs}];
Im[*]:= peakCrossCorrsAft = Table [Max[(catenatedPupilCrossCorrsAft[[n]])[[All, 2]]],
       {n, 1, Length[catenatedPupilCrossCorrsAft]}];
log_{i} = Table[Evaluate@ToExpression[StringJoin["pupModPaired", ToString[usefulR0Is[[n]]]]] = log_{i}
         {usefulROIs[[n]], Flatten@{peakCrossCorrsBef[[n]], peakCrossCorrsAft[[n]]}};,
      {n, 1, Length[usefulROIs]}];
    (***Visualize dFF-pupil correlation data for each ROI for the 2 sessions***)
In[*]:= Manipulate[
     {ListLinePlot[{ToExpression[StringJoin["dffPupCCBef", ToString[n]]], ToExpression[
          StringJoin["dffPupCCAft", ToString[n]]]}, PlotStyle → {Black, Red}, PlotRange → All],
      ToExpression[StringJoin["pupModPaired", ToString[n]]]}, {n, usefulROIs}]
Infel: Export[StringJoin["S:/Imaging/Garrett/FMB208 2PRig/", date, "/", mouse,
       "/PairedAnalysis/", date, "_", mouse, "_pairedROIsPupil.txt"], usefulROIs];
Infer: Table [Export [StringJoin ["S:/Imaging/Garrett/FMB208 2PRig/", date, "/", mouse,
        "/PairedAnalysis/", date, "_", mouse, "_pupilModPaired_ROI", ToString[n], ".txt"],
       ToExpression[StringJoin["pupModPaired", ToString[n]]]], {n, usefulROIs}];
```