

```

    (**Input identifying information**)

In[ ]:= date = ToString[Evaluate[Input["Input the date of the experiment"]]]

In[ ]:= mouse = ToString[Evaluate[Input["Input the mouse identity (e.g. Mouse123)"]]]

In[ ]:= sessionNum = Evaluate[Input["Input the session number"]]

In[ ]:= (**Import the raw pixel intensities from the whisker pad**)
whiskPixIntensities =
    Part[#, 2] & /@ (Drop[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
        mouse, "/Session", ToString[sessionNum], "/", date, "_", mouse, "_", "Session",
        ToString[sessionNum], "_Whiskers/WhiskerPadPixelIntensity/", "Results.csv"]], 1]);

In[ ]:= (**Import the frame times from the camera**)

In[ ]:= frameTimes =
    Drop[Drop[(Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
        "/Session", ToString[sessionNum], "/", date, "_", mouse, "_",
        "Session", ToString[sessionNum], "_CamSync.txt"], "List"]], 16], -1];

In[ ]:= (**Calculate camera frame rate**)
frameRate = Round[Length[frameTimes] / (Last[frameTimes] - First[frameTimes])]

    (**Verify that length of frameTimes = length of whiskPixIntensities, +/- 2**)

In[ ]:= Length[whiskPixIntensities]

In[ ]:= Length[frameTimes]

In[ ]:= (**Match up the length of the cam sync pulses with the length of the image frames**)

In[ ]:= If[Length[frameTimes] > Length[whiskPixIntensities],
    frameTimes = Take[frameTimes, Length[whiskPixIntensities]];
    whiskPixIntensities = Take[whiskPixIntensities, Length[frameTimes]];];

In[ ]:= (**Calculate the motion energy of the whisker pad**)

In[ ]:= whiskMotionEnergy = Abs /@ Differences[whiskPixIntensities];

In[ ]:= whiskMotionEnergyTS = Partition[Riffle[Drop[frameTimes, 1], whiskMotionEnergy], 2];

In[ ]:= whiskMotionInterp = Interpolation[whiskMotionEnergyTS];

In[ ]:= whisker10Hz = Table[whiskMotionInterp[t],
    {t, First[whiskMotionEnergyTS][[1]], Last[whiskMotionEnergyTS][[1]], 0.1}];

In[ ]:= whisker10HzTimeSeries = Partition[Riffle[Table[n, {n, First[whiskMotionEnergyTS][[1]],
    Last[whiskMotionEnergyTS][[1]], 0.1}], whisker10Hz], 2];

In[ ]:= (**Import calcium fluorescence traces**)

In[ ]:= numAdditionalMovies =
    Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/",
        date, "/", mouse, "/Session", ToString[sessionNum], "/Ftraces/"]]] - 2;
    (*Subtract 2 because the first movie was already imported and the
    other file in the directory is an ROI list*)

```

```

In[ ]:= tpFrameTimes =
  Drop[Drop[(Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
    mouse, "/Session", ToString[sessionNum], "/", date, "_", mouse, "_",
    "Session", ToString[sessionNum], "_2PFrameTimes.txt"], "List"]], 16], -1];

In[ ]:= numFramesPerAcq = Round[Length[tpFrameTimes] / (numAdditionalMovies + 1)];

In[ ]:= tpFrameRate = Round[Mean[numFramesPerAcq / Flatten[Differences /@
  ({First[#], Last[#]} & /@ (Partition[tpFrameTimes, numFramesPerAcq]))]]];

In[ ]:= (**For each ROI picked for the session, upload the extracted dF/F0 time series**)

In[ ]:= numROIs =
  Length[FileNames["*", File[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
    mouse, "/Session", ToString[sessionNum], "/dFoverF0TimeSeries/"]]]];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFtimeseries", ToString[n]]] =
  ToExpression /@ Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date,
    "/", mouse, "/Session", ToString[sessionNum], "/dFoverF0TimeSeries/",
    date, "_", mouse, "_Session", ToString[sessionNum], "_",
    "dFoverF0ts_ROI", ToString[n], ".txt"], "List"]; {n, 1, numROIs}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFts", ToString[n]]] =
  Partition[Riffle[(ToExpression[StringJoin["dFFtimeseries", ToString[n]]])][[All, 1]],
    Standardize[(ToExpression[StringJoin["dFFtimeseries", ToString[n]]])][[All, 2]]],
    2]; {n, 1, numROIs}];

In[ ]:= (**Lowpass filter dF/F traces at 10 Hz**)

In[ ]:= Table[
  Evaluate@ToExpression[StringJoin["dFFlpfTS", ToString[n]]] = LowpassFilter[TimeSeries@
    ToExpression[StringJoin["dFFts", ToString[n]]], 10 Hz]; {n, 1, numROIs}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFlpfInt", ToString[n]]] =
  Interpolation[ToExpression[StringJoin["dFFlpfTS", ToString[n]]]]; {n, 1, numROIs}];

In[ ]:= (**Import time periods with no locomotion**)

In[ ]:= quiescentPeriods =
  ToExpression[Import[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/",
    mouse, "/Session", ToString[sessionNum], "/LocomotionData/", date, "_", mouse,
    "_", "Session", ToString[sessionNum], "_QuiescentBouts.txt"], "List"]];

In[ ]:= (**Only keep quiescent periods at least 9 s long**)

In[ ]:= quiescentPeriodsLongBef =
  DeleteCases[Table[If[quiescentPeriods[[n, 2]] - quiescentPeriods[[n, 1]] ≥ 9,
    quiescentPeriods[[n]], Null], {n, 1, Length[quiescentPeriods]}], Null];

In[ ]:= (**For each acquisition bout, determine the quiescent periods
  contained within it and then make sure they're at least 9 s long**)

In[ ]:= acqBoutIntervals =
  Interval /@ ({First[#] + 8, Last[#] - 8} & /@ (Partition[tpFrameTimes, numFramesPerAcq]));

```

```

In[ ]:= testedIntervals = Table[MinMax /@ (IntervalIntersection[acqBoutIntervals[[n]], #] & /@
      (Interval /@ quiescentPeriodsLongBef)), {n, 1, Length[acqBoutIntervals]};

In[ ]:= acceptedQuiescentPeriods = Flatten[
      Table[DeleteCases[Table[If[NumberQ@ (Differences[testedIntervals[[m, n]]][[1]]),
            testedIntervals[[m, n]], Null], {n, 1, Length[testedIntervals[[m]]}],
            Null], {m, 1, Length[testedIntervals]}], 1];

In[ ]:= quiescentPeriodsLong = DeleteCases[
      Table[If[acceptedQuiescentPeriods[[n, 2]] - acceptedQuiescentPeriods[[n, 1]] ≥ 9,
            acceptedQuiescentPeriods[[n]], Null],
            {n, 1, Length[acceptedQuiescentPeriods]}], Null];

In[ ]:= (*****
      *****
      (**In this part, calculate cross-correlograms between dF/F0 and whisker motion**)
      *****
      *****)

In[ ]:= (**Interpolate the 10 Hz sampled whisker time series**)

In[ ]:= whisker10HzInt = Interpolation[whisker10HzTimeSeries];

In[ ]:= (**Compute cross-correlations between whisker and dF/F0 during quiescent
      periods. Compute in a time window of 8 seconds with 100 ms resolution**)

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFwhiskCrossCorr", ToString[roi]]] =
      Mean[DeleteCases[Table[whisker = Table[whisker10HzInt[t],
            {t, quiescentPeriodsLong[[i, 1]], quiescentPeriodsLong[[i, 2]], 0.1}];
            neuron = Table[(ToExpression[StringJoin["dFFlpfInt", ToString[roi]]])[t],
            {t, quiescentPeriodsLong[[i, 1]], quiescentPeriodsLong[[i, 2]], 0.1}];
            If[quiescentPeriodsLong[[i, 2]] - quiescentPeriodsLong[[i, 1]] > 8,
            Join[Reverse[Table[Correlation[whisker, RotateRight[neuron, -n]], {n, 0, 80, 1}]],
            Drop[Table[Correlation[whisker, RotateRight[neuron, n]], {n, 0, 80, 1}], 1]],
            Null], {i, 1, Length[quiescentPeriodsLong]}], Null]], {roi, 1, numROIs}];

In[ ]:= (**For each ROI, make a time series of cross-correlation**)

In[ ]:= crossCorrTimeVals = Table[n, {n, -8, 8, 0.1}];

In[ ]:= Table[Evaluate@ToExpression[StringJoin["dFFwhiskCrossCorrTS", ToString[n]]] =
      Partition[Riffle[crossCorrTimeVals, (ToExpression[
            StringJoin["dFFwhiskCrossCorr", ToString[n]]), 2]], {n, 1, numROIs}];

      (**Visualize whisker-dFF cross correlograms**)

In[ ]:= Manipulate[ListLinePlot[ToExpression[StringJoin["dFFwhiskCrossCorrTS", ToString[roi]]],
      {roi, 1, numROIs, 1}]

In[ ]:= (**Also calculate a summary cross correlation value for each ROI**)

```

```

In[ ]:= Table[Evaluate@ToExpression[StringJoin["summaryDffWhiskCorr", ToString[roi]]] = Mean[
  Table[Correlation[Table[ToExpression[StringJoin["dFFlpfInt", ToString[roi]]][t],
    {t, quiescentPeriods[[n, 1]], quiescentPeriods[[n, 2]], 0.1}], Table[
    whisker10HzInt[t], {t, quiescentPeriods[[n, 1]], quiescentPeriods[[n, 2]], 0.1}]],
    {n, 1, Length[quiescentPeriods]}]]]; {roi, 1, numROIs}];

In[ ]:= (**Export data**)

In[ ]:= Table[Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse, "/Session",
  ToString[sessionNum], "/", "WhiskerData/", date, "_", mouse, "_", "Session",
  ToString[sessionNum], "_dFFwhiskCrossCorr_ROI", ToString[n], ".txt"],
  ToExpression[StringJoin["dFFwhiskCrossCorrTS", ToString[n]]]], {n, 1, numROIs}];

In[ ]:= Table[Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse, "/Session",
  ToString[sessionNum], "/", "WhiskerData/", date, "_", mouse, "_", "Session",
  ToString[sessionNum], "_summaryDffWhiskCorr_ROI", ToString[n], ".txt"],
  ToExpression[StringJoin["summaryDffWhiskCorr", ToString[n]]]], {n, 1, numROIs}];

whiskModROIs = Range[numROIs];

In[ ]:= Export[StringJoin["S:/Imaging/Garrett/FMB208_2PRig/", date, "/", mouse,
  "/Session", ToString[sessionNum], "/", "WhiskerData/", date, "_", mouse, "_",
  "Session", ToString[sessionNum], "_WhiskModROIs", ".txt"], whiskModROIs];

```